# RapidWarm™ Oocyte

## US

## EN: Indication for use

Media for warming of vitrified oocytes.

**Product Description**
SUPPLEMENTED WITH HSA
RapidWarm™ Oocyte contains four solutions for the warming of vitrified oocytes. The solutions consist of a MOPS- buffered medium containing gentamicin as an antibacterial agent, and human serum albumin.
Warm 1™ Oocyte contains sucrose
Warm 2™ Oocyte contains sucrose
Warm 3™ Oocyte contains sucrose
Warm 4™ Oocyte contains no sucrose
For use after warming to +37 °C in ambient atmosphere.

**Storage instructions and stability**
Store in dark at +2 to +8 °C.
RapidWarm™ Oocyte is stable until the expiry date shown on the container and the LOT-specific Certificate of Analysis.
Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the refrigerator. Record opening date on the bottle. Discard excess media no later than two weeks after first opening.

**Directions for use**
The following is the general procedure for using RapidWarm™ Oocyte.
Warning should only be performed by staff trained in warming procedures.
Note: Timing with warming procedures is critical, ensure you follow the protocol precisely.

**Warning:**
Place 1 ml of each of the following media into separate wells of a multi-well plate and warm to 37 °C in ambient atmosphere:

- Warm 1™ Oocyte
- Warm 2™ Oocyte
- Warm 3™ Oocyte
- Warm 4™ Oocyte

NOTE: The recommended volumes should not be changed. Volume changes will affect temperature control in the first warming solution as well as osmolality, which may result in suboptimal oocyte survival.

It is very important to keep the temperature of all warming solutions at 37 °C at all times. If more than one Rapid™ is to be warmed in the same dish, make sure that the temperature reaches 37 °C after warming procedure.

All manipulations of the oocytes are carried out at 37 °C (on the heated stage) in ambient atmosphere. Remove the cryovialde containing the vitrified oocytes from the cryostorage container. Follow the warming instructions for the specific cryovialde used. Immediately, transfer the oocytes into Warm 1™ Oocyte. Important to use 1ml of Warm 1™ Oocyte to ensure the temperature of the medium is not perturbed by successive cryovialde submersions.

Allow the oocytes to fall from the device and sink to the bottom. Leave for 1 minute.

Transfer the oocytes into Warm 2™ Oocyte and let osmolality remain in the solution for 3 minutes. Transfer the oocytes into Warm 3™ Oocyte and let the oocytes remain in the solution for 5 minutes. Transfer the oocytes into Warm 4™ Oocyte and let the oocytes remain in the solution for 5 to 10 minutes.

Following warming, oocytes should be held in an incubator in equilibrated G-1™ PLUS at 37 °C at 6% CO2. The oocytes should be washed several times through this medium.

Perform ICSI 1-2 hours after warming. Continue culture according to standard laboratory practice.

**Specifications**
Sterile filtered SAL 10+
Mouse Embryo Assay (1-cell) [% expanded blastocyst within 96 hours] ≥ 80
Bacterial endotoxins (LAL assay) [EU/mL] < 0.5
LOT specific test results are available on the Certificate of Analysis provided with each delivery.

**Precautions**
Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use.
The long-term safety of vitrification on children born following this method of oocyte cryopreservation procedure has not been established.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife’s IVF media, have not been determined and are uncertain.
The safety and effectiveness of vitrification has not been fully evaluated in human oocytes.

Discard product if bottle integrity is compromised. Do not use RapidWarm™ Oocyte if it appears cloudy. RapidWarm™ Oocyte contains gentamicin and human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBs, HCV, and HTLV I/II RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.
Not for injection.

<b>DESCRIPTION</b>	Sterilized using aseptic processing techniques
<b>TEMPERATURE</b>	Use by - see label
<b>CAUTION:</b>	Consult accompanying documents.
<b>LOT</b>	Catalog number
<b>CE</b>	Batch code
<b>CE</b>	CE Mark (Conformité Européen)

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