

EN: Indication for use

Medium for oocyte retrieval and rinsing (follicle flushing).

Product Description

ASP™ is a bicarbonate and HEPES buffered medium containing heparin and gentamicin as an antibacterial agent.

For use after equilibration at +37°C and 5 % CO₂ atmosphere.

Storage Instructions and stability

Store dark +2 to +8°C.

ASP™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

Equalibrate ASP™ at +37°C and 5 % CO₂.

Rinse the aspiration needle lumen and tubing using ASP™-medium and discard the rinsed fluid. The follicles aspirates may be aspirated individually or collectively. The aspirated follicles must be collected in a sterile, non-tissue culture tube culture.

It is recommended to keep ASP™ at +37°C during the procedure to prevent coagulation of the blood, and to store the collected container to retain pH.

The follicle aspirates should be collected in a sterile, tissue culture disposable vessel and should be examined by the laboratory immediately. If they cannot be examined straight away, the vessels should be tightly sealed and kept at +37°C.

Specifications

Sterile filtered SAL 10³

Mouse Embryo Assay (1-cell) [% expanded blastocyst within 96 hours] ≥ 80

Bacterial endotoxins (LAL assay) [EU/ml] < 0.25

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised.

Do not use ASP™ if it appears cloudy.

ASP™ contains heparin and gentamicin.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not to be used as an injectable product other than follicle flushing during oocytes retrieval.

Caution: Federal (US) law restricts this device to sale only by or on the order of a physician.

Description of ISO Symbols

Sterilized using aseptic processing techniques

Temperature limit

Do not reuse, discard after procedure

Use by - see label

Caution: Consult accompanying documents.

Catalog number

Batch code

CE Mark (Conformité Européenne)

CHINA**CN: 取卵液 ASP™ 产品说明书****适应证**

本培养液用于取卵和冲洗（卵泡）过程。

产品型号及规格

型号：见标签

规格：见标签

产品注册证及技术要求编号

产品注册证编号：见标签

产品批号要求说明：见标签

产品描述

ASP™是含有肝素和抗凝剂庆大霉素的碳酸氢盐和乙酰胺（HEPES）缓冲的培养液。

在+37°C和CO₂浓度为5%的大气环境中平衡后使用。

培养液必须在开封后应再继续储藏。使用后丢弃并密封保存。

使用方法

在+37°C和CO₂浓度为5%的大气环境下平衡。

使用ASP™冲取针腔针孔和连接管，丢弃冲洗后的培养液。卵泡可单个或成簇收集。取出的卵泡必须在无菌的培养液中平衡后使用。

建议每个操作过程中使用恒温以保持ASP™培养液。

卵泡液必须在平衡后应再继续储藏。使用后丢弃并密封保存。

注意事项

如果出现瓶子的整体完整性受损，请丢弃该产品。如果发现产品包装盒未用尽，请丢弃该产品。

ASP™含有肝素和庆大霉素。

重复使用可能致微生物污染或产品性能改变。

Vitrolife建议在使用本产品时，严格执行无菌操作技术。

关于使用IVF培养液产品，包括Vitrolife的IVF培养液产品，在生殖毒性和发育毒性上的风险还没有论断。

本产品为取卵期间的冲洗所用，不得作为注射用制剂。

注意：产品必须便符合医疗部门相关操作规范和法规，仅限于治疗的医务人员使用。

ISO符号描述

经无菌处理

温度限制

切勿再次使用

使用期限/见标签

注意：参考随附文件

型号

批号

CE 标志 (Conformité Européenne)

EUROPE**EN: Indication for use**

Medium for oocyte retrieval and rinsing (follicle flushing).

Product Description

ASP™ is a bicarbonate and HEPES buffered medium containing heparin and gentamicin as an antibacterial agent.

For use after equilibration at +37°C and 5 % CO₂ atmosphere.

Storage Instructions and stability

Store dark +2 to +8°C.

ASP™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

Equalibrate ASP™ at +37°C and 5 % CO₂.

Rinse the aspiration needle lumen and tubing using ASP™-medium and discard the rinsed fluid.

The follicles aspirates may be aspirated individually or collectively. The aspirated follicles must be collected in a sterile, non-tissue culture tube culture tube.

It is recommended to keep ASP™ at +37°C during the procedure to prevent coagulation of the blood, and to store the collected container to retain pH.

The follicle aspirates should be collected in a sterile, tissue culture disposable vessel and should be examined by the laboratory immediately. If they cannot be examined straight away, the vessels should be tightly sealed and kept at +37°C.

Specifications

Sterile filtered SAL 10³

Mouse Embryo Assay (1-cell) [% expanded blastocyst within 96 hours] ≥ 80

Bacterial endotoxins (LAL assay) [EU/ml] < 0.25

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised.

Do not use ASP™ if it appears cloudy.

ASP™ contains heparin and gentamicin.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not to be used as an injectable product other than follicle flushing during oocytes retrieval.

Caution: Federal (US) law restricts this device to sale only by or on the order of a physician.

Description of ISO Symbols

Sterilized using aseptic processing techniques

Temperature limit

Do not reuse, discard after procedure

Use by - see label

Caution: Consult accompanying documents.

Catalog number

Batch code

CE Mark (Conformité Européenne)

ASPs™ 对于卵泡液的贮存和稳定性

ASPs™ 是稳定的，直到其过期日期。

ASPs™ 在 +37°C 和 5% CO₂ 下稳定。

ASPs™ 不应存储在过期的日期。

ASPs™ 不应存储在过期的