HYASE™-10X

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US

EN: Indication for use

Medium for removal of cumulus cells.

Product Description

HYASE™ is a physiological salt buffer containing hyaluronidase, human serum albumin and penicillin G as an antibacterial agent.

For use after dilution 1:10 with either of the following two alternatives:

1) G-MOPSTM PLUS or G-MOPSTM supplemented with G-MM™ or HSA-solution™ and equilibration at +37°C and ambient atmosphere.

2) G-GAMETETM and equilibration at +37°C and 6 % CO₂ atmosphere.

The concentration of hyaluronidase in HYASE™ is 800 IU/ml before dilution.

Storage instructions and stability

Store dark at +2 to +8 °C.

HYASE™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

Dilute HYASE™ either with supplemented G-MOPSTM/G-MOPSTM PLUS/G-GAMETETM to a ratio of 1:10, mix well.

Prepare and label sterile, non toxic multi-well dishes and prepare 3 wells of either supplemented G-MOPSTM/G-MOPSTM PLUS/G-GAMETETM for every well of HYASE™.

Pre-warm the multi well dishes with diluted HYASE™ and supplemented G-MOPSTM/G-MOPSTM PLUS for 15 minutes.

If G-GAMETETM is used as diluting medium, equilibrate the dishes at +37 °C and 6 % CO₂ for 3 hours.

Using a sterile, large bore pipette, place 3–5 oocytes into the HYASE™. Gently pipette the HYASE™ and the oocytes to disperse the cumulus cells. Do not expose oocytes to HYASE™ solution for more than 30 seconds.

Use a sterile pipette to transfer oocytes from HYASE™ to the first washing well. Using a fine bore denudation pipette, aspirate each oocyte singly up and down to remove the corona. The diameter of the pipette should be slightly larger than that of the oocyte. Using a pipette that is too narrow or has a jagged edge may damage the oocyte.

Use the remaining two washing wells to rinse the oocyte. Transfer oocyte to new supplemented G-MOPSTM/G-MOPSTM PLUS/G-GAMETETM or other equivalent droplets for oocyte handling or, if needed, to culture medium for further incubation and maturation.

Specifications

Bacterial endotoxins (LAL assay) [EU/ml]	SAL 10 ³	< 0.5*
*after diluting 1:10 in G-MOPSTM PLUS		
Cumulus Cell Removal Assay	PASS	
Mouse Embryo Assay	Not performed	
LOT specific test results are available on the Certificate of Analysis provided with each delivery.		

Precautions

Discard product if bottle integrity is compromised. Do not use HYASE™ if it appears cloudy.

HYASE™ contains hyaluronidase, human serum albumin and penicillin G.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBs, HCV, and HTLV III and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.
Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

Description of ISO Symbols

 STERILIZABLE	STERILIZED using aseptic processing techniques
 Temperature limit	Temperature limit
 Do not re-use, discard after procedure.	Do not re-use, discard after procedure.
 Use by - see label.	Use by - see label.
 Caution: Consult accompanying documents.	Caution: Consult accompanying documents.
 Catalog number	Catalog number
 LOT	LOT
 CE Mark (Conformité Européen)	CE Mark (Conformité Européen)

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