RapidWarm™ Blast

US Vitrolife

US

EN: Indication for use

RapidWarm™ Blast is indicated for warming of vitrified human blastocyst stage embryos.

Product Description

Media for warming of vitrified human blastocyst stage embryos. The media consists of a MOPS buffered solution containing gentamicin as an antibacterial agent and human serum albumin.

SUPPLEMENTED WITH HSA

RapidWarm™ Blast contains three media for the warming of vitrified human blastocyst stage embryos. The media consists of a MOPS buffered solution containing gentamicin as an antibacterial agent and human serum albumin. Warm 1™ Blast contains sucrose as cryoprotectant.

Warm 2™ Blast contains sucrose as cryoprotectant. Warm 3™ Blast contains no cryoprotectants. For use after warming to +37 °C.

Storage instructions and stability

Store dark at +2 to +8 °C.

RapidWarm™ Blast is stable until the expiry date shown on the bottle labels and the LOT specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of procedure.

Use of storage device

Use a legally marked storage device indicated for use in human blastocyst vitrification procedures, in order to assure sufficient cooling and warming rates.

Use a closed storage system to prevent the potential risk of viral or other contamination of samples.

Perform the actual vitrification and warming according to the instructions for Use for the storage device.

Directions for use

Caution: The user should read and understand the Directions for Use, Warnings and Precautions, and be trained in the correct procedure before using Vitrolife RapidVi™ Blast and RapidWarm™ Blast kits for vitrification and warming of human blastocyst stage embryos.

The following is the general procedure for using RapidWarm™ Blast.

Note that the timing has to be accurate.

Warning: Place 0.5-1 ml of each of the following media into separate wells of a 4-well plate and warm to 37 °C.

• Warm 1™ Blast

• Warm 2™ Blast

• Warm 1™ Blast

All manipulations of the blastocysts are carried out at 37 °C (on the heated stage).

Remove the cryodevice containing the vitrified blastocysts from the cryostorage container. Follow the warming instructions for the specific cryodevice used.

Immediately after warming, transfer blastocysts into Warm 1™ Blast.

Allow the blastocysts to fall from the device and sink to the bottom. Leave for 2 min.

Transfer the blastocysts into Warm 2™ Blast and let the blastocysts remain in the solution for 3 min.

Transfer the blastocysts into Warm 3™ Blast and let the blastocysts remain in the solution for 5-10 min.

Rinse the blastocysts in culture media several times and continue culture according to laboratory protocol. Embryo(s) can be cultured using G-2™ or EmbryoGlu® for this procedure.

Specifications

Sterile filtered SAL 10⁺ Mouse Embryo Assay (1-cell) [% re-expanded blastocyst 24 hours post-test] ≥ 70

Bacterial endotoxins (LA assay) [EU/mL] < 0.5

LOT-specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use.

The long-term safety of vitrification and/or blastocyst collapse on children born following this method of embryo cryopreservation procedure has not been established.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

The safety and effectiveness of vitrification has not been fully evaluated in human embryos that have not yet reached the blastocyst stage of development.

Discard product if bottle integrity is compromised. Do not use RapidWarm™ Blast if it appears cloudy.

RapidWarm™ Blast contains gentamicin and human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBs, HCV and HTLV III and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis.

RapidWarm™ Blast media can offer assurance that products derived from human blood will not transmit infectious agents.

Re-use may result in microbiological contamination and/or property changes in the media.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

Not for injection.

Description of ISO Symbols

STERILE	Sterilized using aseptic processing techniques	
—④—	Temperature limit	
—⑤—	Do not re-use, discard after procedure.	
—⑥—	Use by - see label.	
—⑦—	Caution: Consult accompanying documents.	
—⑧—	Catalog number	
—⑨—	Batch code	
CE 200	CE Mark (Conformité Européen)	

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