SpermFreeze Solution™

Vitrolife

US

EN: Indication for use

Medium for cryopreservation of human sperm.

Product Description

SpermFreeze Solution™ is a bicarbonate and MOPS-based medium containing gentamicin as an antibiotic agent and human serum albumin.

For use after equilibration to +20 °C ± 5°C.

Storage instructions and stability

Store dark at +20 °C.

SpermFreeze Solution™ is stable until the expiry date shown on the container labels and the LOT specific Certificate of Analysis.

Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the bottle. Record dating date on the bottle. Discard excess media no later than two weeks after first opening.

Directions for use

The following is the general procedure for using SpermFreeze Solution™. If you have not used this product before, consult the cryopreservation method is essential, to ascertain sperm survival.

When using straws:

Frosting

Allow the semen to liquefy for at least 10 minutes prior to preparation and ensure that both semen and SpermFreeze Solution™ are at room temperature.

Measure the total volume of semen and carry out semen analysis as required.

Dilute with equal volume of semen and SpermFreeze Solution™. Semen and SpermFreeze should be added slowly to the semen to ensure the mixture carefully tilted after each drop added. Close the lid tightly and turn the tube upside down 20 times, being careful not to create bubbles.

The semen mixture is left at room temperature for 10 minutes.

Mark 5 mL straws with patient ID and load the semen mixture into straws. Ensure that some air is left in the lower part of the straw for sealing as well as to allow expansion of the semen mixture during freezing. Seal the straw horizontally on a 1-cm strip of paper.

Place the straws horizontally on a 1-cm strip of paper board in a liquid nitrogen bath to allow for freezing in vapour phase. Leave for 30 minutes*.

Transfer the straws quickly into liquid nitrogen and store at -196 °C.

* Optimal results can be performed using a slow-freeze machine programmed for sperm freezing.

Frosting program for machine:
Start temperature: +20 °C
1...5 °C/min to -8 °C
2. Hold at -8 °C
3...10 °C/min to -25 °C
4...25 °C/min to -150 °C

Thawing

Work aseptically.

Remove straws from -196 °C and place them in a water bath at +20 °C for 10 minutes.

Wipe the straws dry with a clean paper towel and open according to instructions from the manufacturer of the straws.

Expel semen mixture into these test tubes and dilute with equal amount of equilibrated G-IV™ PLUS. G-IV™ PLUS should be added dropwise to the semen and the solution carefully mixed after each addition.

Continue with gradient separation according to the G-series manual.

When using cryovials:

Frosting

Allow the semen to liquefy for at least 10 minutes prior to preparation and ensure that both semen and SpermFreeze Solution™ are at room temperature.

Measure the total volume of semen and carry out semen analysis as required.

Use with equal volume of semen and SpermFreeze Solution™. Semen and SpermFreeze should be added slowly to the semen and the mixture carefully tilted after each drop added. Close the lid tightly and turn the tube upside down 20 times, being careful not to create bubbles.

The semen mixture is left at room temperature for 10 minutes.

Mark 5 mL cryovials with patient ID and load the semen mixture into cryovials. Do not fill cryovials completely to allow for expansion of the semen mixture during freezing.

Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the bottle. Record dating date on the bottle. Discard excess media no later than two weeks after first opening.

The following is the general procedure for using a slow-freeze machine programmed for sperm freezing.

Frosting program for machine:
Start temperature: +20 °C
1...5 °C/min to -8 °C
2. Hold at -8 °C
3...10 °C/min to -25 °C
4...25 °C/min to -150 °C

Thawing

Work aseptically.

Remove cryovials from -196 °C and place them in a water bath at +20 °C for 10 min.

Wipe the cryovials dry with a clean paper towel.

Transfer the semen mixture to clean test tubes.

Dilute the semen mixture with equal amount of equilibrated G-IV™ PLUS. G-IV™ PLUS should be added dropwise to the semen and the mixture carefully mixed after each addition.

Continue with gradient separation according to the G-series manual.

Specifications

STERILE filtered

HUMAN SERUM SURVIVAL ASSAY [%]

Bacterial endotoxins (LAL assay) [EU/ml]

Mouse Embryo Assay not performed

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised. Do not use SpermFreeze Solution™ if it appears cloudy.

SpermFreeze Solution™ contains human serum albumin and gentamicin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBc, HBe, and HCV and non-reactive for hepatitis C antibody. No test methods are available for HIV-1 and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's VVF media, have not been determined and are uncertain.

No injection.

Description of ISO Symbols

STERILE A Sterilized using aseptic processing techniques

-80°C Temperature limit

Use by - see label.

Caution Consult accompanying documents.

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