

Cell Culture

Hypure™ Oil

- Light and Heavy, 2 types are available.
- High-purity paraffin oil without RNA detection.
- Prepared oil which is washed twice with ultrapure water.
- Light shielding glass is used to prevent deterioration of quality by light. Air tightness is kept in the bottle container.



Enlargement

	Order Number	Code	Contents
Hypure Oil Heavy	93621	MOHV-100	100mL
Hypure Oil Light	93521	MOLH-100	

COMPONENTS

High-purity Paraffin Oil

Specific Gravity	Light	Heavy
	0.838g/mL, 25°C	0.862g/mL, 25°C

QUALITY CONTROL

Sterile Filtration / $POV \leq 0.1$ meq/kg / Readily Carbonizable Substances / Endotoxin < 0.25 EU/mL / Mouse Embryo Assay $\geq 80\%$ / Sterility test

Storage: 8-30°C / Shelf Life: 12 months